

<b>College/Unit:</b>	<b>Department of Mechanical Engineering</b>
<b>Procedure Title:</b>	<b>Standard Operating Procedure #TE0003 Decontamination Procedures</b>

<b>Author:</b>		
<b>Departmental Assistant</b>	- <b>Doug Bitner</b>	
<b>Ph.D. Candidate</b>	- <b>Liqun Ning</b>	
Print Title	Print Name	Signature

<b>Authorized By:</b>		
<b>Professor</b>	<b>Daniel Chen</b>	
Print Title	Print Name	Signature

## Contents

Revision History .....	1
1 Purpose .....	2
2 Applicable To .....	2
3 Scope .....	3
4 Definitions .....	3
5 Training Required .....	3
6 Safety .....	3
7 Equipment and Materials .....	3
8 Procedure .....	4
9 References .....	5

## Revision History

Revisions to this procedure are to be documented in Table 1, Revision History.

Table 1: Revision History

<b>Document Section</b>	<b>Details of Amendments</b>	<b>Date</b>	<b>Author (Initials)</b>



### **3 Scope**

This procedure covers equipment and materials in Room 1A26.

### **4 Definitions**

SOP: Standard Operating Procedure

PPE: Personal Protective Equipment

SDS: Safety Data Sheets

### **5 Training Required**

Laboratory Safety

WHMIS2015

Biosafety

### **6 Safety**

Check SDS for safe handling procedures of the sterilizing and decontamination solutions used.

PPE requirements

- Safety glasses
- Lab Coat
- Appropriate gloves for the sterilizing and decontamination solutions used.

### **7 Equipment and Materials**

The equipment to be decontaminated includes the 3D-Bioplotter, clean benches, incubator, scaffolding material and cultured cells. The chemicals used to decontaminate or sterilize specific equipment or materials are outlined in the following procedure.

## 8 Procedure

This section provides the steps involved in carrying out the decontamination and sterilization procedures.

### 1) 3D-Bioplotter and Terra Universal Clean Bench sterilization process

- a. Turn on both clean bench fans;
- b. Using a tissue and a solution of 70% ethanol, wipe and clean the surfaces of the Bioplotter, clean bench and work area. During this procedure, make sure that the Bioplotter is turned off;
- c. After 10 minutes, repeat step b;
- d. Keep fans on for at least 1 hour, then the Bioplotter work can be started.

### 2) Envirco Clean Bench sterilization process

- a. Turn on the clean bench fan;
- b. Using a tissue and a solution of 70% ethanol, wipe and clean the surfaces of the clean bench and work area.
- c. After 10 minutes, repeat step b;
- d. Keep fans on for at least 1 hour, then work can be started.

### 3) Tissue culture incubator decontamination

- a. Turn off the incubator and the CO<sub>2</sub> tank valve and take out all of the shelves;
- b. Wipe the interior of the incubator completely using a tissue and a solution of 50% bleach;
- c. Wipe the incubator shelves completely using a tissue and a solution of 50% bleach;
- d. Wipe the interior of the incubator using a tissue and deionized water;
- e. Wipe the shelves using a tissue and deionized water;
- f. Wipe the interior of the incubator completely using a tissue and solution of 70% ethanol;
- g. Wipe the shelves using a tissue and a solution of 70% ethanol;
- h. Put the shelves back into the incubator, turn the regulator of the CO<sub>2</sub> tank to the correct pressure, and turn on the incubator.

### 4) Material sterilization

- a. Use a 0.22 µm bottle-top filter to filter the material solution;
- b. Freeze and then freeze-dry the filtered material solution;
- c. Dissolve freeze-dried material in sterilized solvent in a sterilized container;
- d. The material should now be ready for use.

## **5) Decontamination of cells**

- a. Once cultured cells are contaminated, cells cannot be used anymore.
- b. Add 10% bleach solution to cell culture medium to kill both bacteria and cells.
- c. Dispose of mixture down the sink and decontaminate the sink.

## **9 References**

Canadian Biosafety Handbook, Second Edition, May 26, 2016